



Product Information

LA974

Micropipette Tips Capacity: 200 µL

Product Number

Packing

LA974 : 1 X 1000 No.

Product Information

Product Code : LA974
Product Name : Micropipette Tips Capacity: 200 µL

Technical Specification

Description : Autoclavable transparent nuclease free, pyrogen free micropipette tips.
Clarity : Clear transparent without any particles and Fibers.
Inner Dimension : Approx 5.3 mm
Capacity : 200 µL
Product alignment : Straight without any bending.
Compatible Pipette Code : LA599, LA600, LA604, LA612, LA613, LA615, LA618, LA959, LA966.
Temperature tolerance test : Withstands steam sterilization at 121°C for 30 minutes and not to be used in dry heat sterilization.
Use : Used to dispense small amount of liquid.

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods