



# Technical Data

## LA975 1000 µL Micropipette Tips (Transparent)

Product Number	Packing
LA975	: 1 X 1000 No.

### Product Information

Product Code	: LA975
Product Name	: 1000 µL Micropipette Tips (Transparent)

### Technical Specification

Description	: Autoclavable transparent nuclease free, pyrogen free micropipette tips.
Clarity	: Clear transparent without any particles and Fibers.
Inner Dimension	: Approx 7.4 mm
Capacity	: 1000 µL
Product alignment	: Straight without any bending.
Compatible Pipette Code	: LA606, LA607, LA608, LA614, LA616, LA960, LA963.
Temperature tolerance test	: Withstands steam sterilization at 121°C for 30 minutes and not to be used in dry heat sterilization.
Use	: Used to dispense small amount of liquid.

### Transport Information

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods

#### Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

HiMedia Laboratories Pvt. Ltd. A-516, Swastik Disha Business Park, Via Vadhani Ind. Est., LBS Marg, Mumbai-400086, India.  
Customer care No.: 022-6147 1919 Email: [techhelp@himedialabs.com](mailto:techhelp@himedialabs.com) Website: [www.himedialabs.com](http://www.himedialabs.com)