



HiCrome™ Cronobacter Isolation Agar(CCI Agar)

M2062I

Intended Use

Recommended for the isolation and identification of *Cronobacter sakazakii* from food products. The composition and performance of this media are as per specifications laid down in ISO /TS 22964: 2017. It can also be used for clinical samples.

Composition**

Ingredients	Gms / Litre
Tryptone#	7.000
Yeast extract	3.000
Sodium chloride	5.000
Sodium deoxycholate	0.250
5-Bromo-4-chloro-3-indolyl α -D-glucopyranoside	1.500
Ammonium iron(III) citrate	1.000
Sodium thiosulfate	1.000
Agar	15.000
Final pH after sterilization (at 25°C)	7.3 \pm 0.2

**Formula adjusted, standardized to suit performance parameters

- Equivalent to Tryptic digest of casein

Directions

Suspend 32.4 grams in 1000 ml purified / distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Enterobacter species are widely distributed in nature occurring in fresh water, soil, sewage, plants, vegetables, animal and human faeces. **Cronobacter sakazakii* has been closely associated with neonatal meningitis and sepsis (1). HiCrome™ Cronobacter isolation Agar is recommended by ISO Committee for the isolation and identification of **C.sakazakii* from food samples (2). The chromogenic substrate (5-Bromo-4-chloro-3-indolyl α -D-glucopyranoside) is cleaved specifically by **C.sakazakii* resulting in the formation of blue green colonies. Other organisms, which do not cleave this substrate, produce colourless colonies. Tryptone and yeast extract provides nitrogenous and carbonaceous compounds, long chain amino acids, vitamins and other essential growth nutrients. Sodium chloride helps in maintaining the osmotic equilibrium of the medium. Sodium deoxycholate inhibits the accompanying gram-positive flora.

Key: *: Formerly known as *Enterobacter sakazakii*

Type of specimen

Clinical samples- blood, urine , cerebrospinal fluid ; Food samples.

Specimen Collection and Handling

For clinical samples, follow appropriate techniques for sample collection and processing as per guidelines (2,3).

For food samples, follow appropriate techniques for sample collection and processing as per guidelines (5).

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

In Vitro diagnostic Use. Read the label before opening the container. Wear protective gloves /protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

1. Slight variation in colour may be observed depending on enzyme production by organism and substrate utilization from the medium.
2. Some species may show poor growth due to nutritional variations.
3. Further biochemical tests must be carried out for confirmation.

Please refer disclaimer Overleaf.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored recommended temperature.

Quality Control

Appearance

Cream to yellow to pink homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Yellow coloured, clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 3.24% w/v aqueous solution at 25°C. pH : 7.3±0.2

pH

7.10-7.50

Cultural Response

Cultural characteristics observed after an incubation at 41.5±1°C for 24±2 hours.

Organism	Inoculum (CFU)	Growth	Recovery	Colour of Colony
<i>Cronobacter sakazakii</i> ATCC 29544 (00214*)	50-100	good-luxuriant	≥50%	blue-green
<i>Cronobacter mytjensii</i> ATCC 51329 (00213*)	50-100	good-luxuriant	≥50%	blue-green
<i>Enterobacter cloacae</i> ATCC 13047 (00083*)	50-100	good-luxuriant	≥50%	colourless without green or blue green colour
<i>Staphylococcus aureus</i> <i>subsp. aureus</i> ATCC 25923 (00034*)	≥10 ⁴	inhibited	0%	
<i>Staphylococcus aureus</i> <i>subsp. aureus</i> ATCC 6538 (00032*)	≥10 ⁴	inhibited	0%	

Key: (*) Corresponding WDCM numbers

Storage and Shelf Life

Store between 15-25°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (2,3).

Reference

- International Organization for Standardization. Microbiology of the food chain- Horizontal method for the detection of *Cronobacter* spp. Draft ISO/ TS 22964, 20176 (E).
- Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
- Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
- Muytjens H. L., Zanen H. C., Sonderkamp H. J. et al, J. Clin Microbiol 18:115-120, 1983.
- Salfinger Y., and Tortorello M.L., 2015, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.

	In vitro diagnostic medical device
	CE Marking
	Storage temperature
	Do not use if package is damaged
	HiMedia Laboratories Pvt. Limited, 23 Vadhani Industrial Estate, LBS Marg,Mumbai-86,MS,India
	CE Partner 4U ,Esdoornlaan 13, 3951 DB Maarn The Netherlands, www.cepartner4u.eu

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.