



Minimal Broth, Davis

M389

Minimal Broth, Davis is recommended for the isolation and characterization of nutritional mutants of *Escherichia coli*.

Composition**

Ingredients	Gms / Litre
Dextrose	1.000
Ammonium sulphate	1.000
Dipotassium phosphate	7.000
Monopotassium phosphate	2.000
Sodium citrate	0.500
Magnesium sulphate	0.100
Final pH (at 25°C)	7.0±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 11.6 grams in 1000 ml distilled water. Heat if necessary to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and dispense as desired.

Principle And Interpretation

Nutritional mutants of *Escherichia coli* obtained by the exposure of wild type *E. coli* to ultra violet light need a nutritionally complete medium to grow. Minimal media can be supplemented with the desired additives to study nutritional characters of the nutritional mutants. Minimal media are the formulations of Davis (1) as described by Lederberg (2). Minimal media contain the necessary nutrients only for the growth of wild type *E. coli* strains. By the random isolation method described by Lederberg, nutritional mutants derived from irradiated cultures of wild type *E. coli* can be isolated (2). These mutants can also be isolated by the use of penicillin as described by Davis and Lederberg (1). *Bacillus subtilis* mutants can be isolated by these techniques and by the penicillin technique also, as described by Nester et al (3).

Dextrose is an energy source. Dipotassium and monopotassium phosphates provide buffering to the medium. Magnesium sulphate and ammonium sulphate are sources of ions that simulate metabolism.

The nutritional supplements to be added to minimal medium depend upon the type of mutant to be screened as for amino acids, vitamins, nucleic acids or other substances. This can be achieved by addition of vitamin assay casamino acids plus tryptophan or a mixture of water soluble vitamins, yeast or nucleic acid extracts.

A cell suspension of wild type *E. coli* is irradiated and cultured on Minimal Agar (M512) supplement with all the necessary growth requirements. This will allow growth of both wild type cells (prototrophs) and mutant cells. The selected colonies are then added to Minimal Broth, Davis and Minimal Broth Davis supplemented with the growth requirements and incubated at 35°C for 24 hours. Growth in the Minimal Broth supplemented with growth requirements and no growth in Minimal Broth indicates a mutant for that particular factor.

Quality Control

Appearance

White to cream homogeneous free flowing powder

Colour and Clarity of prepared medium

Colourless clear solution in tubes

Reaction

Reaction of 1.16% w/v aqueous solution at 25°C. pH : 7.0±0.2

pH

6.80-7.20

Cultural Response

M389: Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours .

Organism	Inoculum (CFU)	Growth
<i>Escherichia coli</i> ATCC 13762	50-100	luxuriant
<i>Escherichia coli</i> ATCC 23724	50-100	luxuriant

Storage and Shelf Life

Store below 30°C in tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Reference

1. Davis B. D., 1949, Proc. Natl Acad. Sci, 35:1.
2. Lederberg J., 1950, Methods in Med. Res., 3:5.
3. Nester E. W., Schafer M. and Lederberg J., 1963, Genetics, 48:529.

Revision : 2 / 2015

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.