

MB056 2-Butanol, For Molecular Biology**Product Number Packing**

MB056 : 500ML

Product Information

Product Code : MB056
Product Name : 2-Butanol, For Molecular Biology
Molecular Formula : C₄H₁₀O
Molecular Weight : 74.12
CAS No. : 78-92-2
EC No. : 201-158-5
HS Code : 2905 14 90
EC Index No. : 603-127-00-5
Shelf Life : 4 years

Technical Specification

Appearance : Colourless clear liquid.
Solubility : Miscible with water.
DNases : None detected
RNases : None detected
Assay : min. 99.00 %

GHS Safety Information

Hazard Statement(s) : H226-H319-H335-H336
Precautionary Statement(s) : P261-P305+P351+P338
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 10-36/37-67
S-Phrase(s) : 13-24/25-26-46-7/9
WGK : 1
Flash Point(°F) : 80.6 °F
Flash Point(°C) : 27 °C
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

UN No. : 1120



Product Information

Class	:	3
Packaging Group	:	3
Marine Pollutant	:	No
ADR/RID	:	1120 3/PG 3
IMDG	:	1120 3/PG 3
IATA	:	1120 3/PG 3

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein