

MB121 3,3',4,4'-Tetraaminobiphenyl tetrahydrochloride, For Molecular Biology

Product Number	Packing
MB121	: 1G
MB121	: 5G

Product Information

Product Code	: MB121
Product Name	: 3,3',4,4'-Tetraaminobiphenyl tetrahydrochloride, For Molecular Biology
Synonym	: 3,3',4,4'-Biphenyltetramine tetrahydrochloride; DAB
Molecular Formula	: C ₁₂ H ₁₄ N ₄ .4HCl
Molecular Weight	: 360.11
CAS No.	: 7411-49-6
EC No.	: 231-018-9
HS Code	: 2921 59 90
Shelf Life	: 4 years

Technical Specification

Appearance	: Grey to tan to purple-brown powder
Solubility	: 33.3 mg soluble in 1 mL of water
DNases	: None detected
RNases	: None detected
FTIR	: Matches with the standard pattern
Assay (AT)	: 97.00 - 102.00%

GHS Safety Information

Hazard Statement(s)	: H302-H319-H341-H350
Precautionary Statement(s)	: P201-P280-P301+P312+P330-P305+P351+P338-P308+P313-P337+P313
Signal Word	: Danger
Hazard Pictogram(s)	



Risk and Safety Information

R-Phrase(s)	: 45-68
S-Phrase(s)	: 53-45
WGK	: 3
RTECS	: DV8753000
Storage Temperature(°C)	: Store at 2 - 8°C
Hazard Symbol(s)	



Transport Information

Marine Pollutant	: No
------------------	------



Product Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein