

## **MB145 DEAE Dextran, For Molecular Biology**

### **Product Number Packing**

---

MB145 : 100G  
MB145 : 25G

### **Product Information**

---

Product Code : MB145  
Product Name : DEAE Dextran, For Molecular Biology  
CAS NO. : 9015-73-0  
Average Molecular Weight : 5,00,000  
Shelf Life : 4 years

### **Technical Specification**

---

Appearance : White to cream hygroscopic powder  
Solubility : 100 mg soluble in 1 mL of water  
pH (5% in water at 25°C) : 4.00 - 7.00  
DNases : None detected  
RNases : None detected  
Loss on drying (at 105°C, 2hr) : <= 7.00%  
Nitrogen content : 3.00 - 3.50%

### **Risk and Safety Information**

---

WGK : 3  
Storage Temperature(°C) : Store below 30°C

### **Transport Information**

---

Marine Pollutant : -  
ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods