

## MBT106-20µg

## Campylobacter jejuni Genomic DNA


HiMedia's *Campylobacter jejuni* Genomic DNA is an ultrapure DNA extracted from *Campylobacter jejuni*. It is used as positive control in *in vitro* diagnostics techniques based on nucleic acid amplification. The vial contains upto 20µg of the purified DNA.

Upon receipt, the product should be stored at -70°C for long-term storage. When using the product, it is recommended that the DNA should be stored in aliquots at -20°C. The product is stable for upto 30 months.

### Quality Control:

Functionally tested in amplification.

Lane M 1 2



Lane 1: DNA purified from *Campylobacter jejuni*

Lane 2: 16S rRNA PCR product of the purified *Campylobacter jejuni* DNA

### Technical Assistance

At HiMedia's, we pride ourselves on the quality and availability of our technical support. For any kind of technical assistance, mail at [mb@himedialabs.com](mailto:mb@himedialabs.com).

PIMBT106\_O/0514

MBT106-00

### Disclaimer :

User must ensure of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

CE A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai - 400 086, India Customer Care No.: 022-6147 1919 Email: [techhelp@himedialabs.com](mailto:techhelp@himedialabs.com)



#### Registered Office :

23, Vadhani Industrial Estate, LBS Marg,  
Mumbai - 400 086, India.  
Tel. : (022) 4017 9797 / 2500 1607  
Fax : (022) 2500 2286

#### Commercial Office

A-516, Swastik Disha Business Park,  
Via Vadhani Indl. Est., LBS Marg,  
Mumbai - 400 086, India

Tel: 00-91-22-6147 1919  
Fax: 6147 1920, 2500 5764  
Email : [info@himedialabs.com](mailto:info@himedialabs.com)  
Web : [www.himedialabs.com](http://www.himedialabs.com)

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein