
Banana powder

PCT0402

Plant Culture Tested

Product Identifier

HS Code : 1106 30 10
Storage : Below 30°C
Shelf life : 2 years

Technical Specification

Appearance : White to tan homogeneous free flowing powder
Solubility : Insoluble in water
pH (3% in water at 25°C) : 4.00 - 6.00
Water (K.F.) : <= 8.00%
Sulphated ash : <= 3.00%
Plant Tissue Culture Test : Passes

Safety Information

UN No. : Not dangerous goods
Class : -
Packing Group : -
WGK : 3