

PCT0414 Coconut water, Sterile filtered, Plant Culture Tested

Product Number

Packing

PCT0414	: 1X100ML
PCT0414	: 5X100ML

Product Information

Product Code	: PCT0414
Product Name	: Coconut water, Sterile filtered, Plant Culture Tested
Shelf Life	: 7 months

Technical Specification

Appearance	: Colourless to slight yellow tint liquid
Clarity	: Clear solution
Solubility	: Miscible with water
pH (as such)	: 4.00 - 6.00
Culture response	: Cultures conditions - Incubation period (5wks), Relative humidity (60±2%), Temperature (25±2°C), Photoperiod Day: Night in hours (16:8)
Shoot culture	: No structural deformity observed, actively growing shoots, no toxicity to shoots
Sterility	: No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

Risk and Safety Information

Storage Temperature(°C)	: Store at -20°C
-------------------------	------------------

Transport Information

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods