
Thiabendazole, Plant Culture Tested

PCT0819


Product Identifier

CAS No.	:	148-79-8
EC No.	:	205-725-8
Molecular Formula	:	C ₁₀ H ₇ N ₃ S
Molecular Weight	:	201.25
HS Code	:	2934 10 00
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to yellow powder
Solubility	:	33.3 mg soluble in 1 mL of dimethylformamide
FTIR	:	Matches with the standard pattern
Melting range	:	298 - 301°C
Assay (NT/GC)	:	min. 98.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H410
Precautionary Statement(s)	:	P273- P501
UN No.	:	3077
Class	:	9
Packing Group	:	III
RTECS	:	DE0700000
WGK	:	2