

**PCT1115 Amoxicillin:clavulanic acid potassium salt (5:1) (Augmentin),
Plant Culture Tested**

Product Number

Packing

PCT1115 : 2G

Product Information

Product Code : PCT1115
Product Name : Amoxicillin:clavulanic acid potassium salt (5:1) (Augmentin), Plant Culture Tested
Synonym : Augmentin
CAS No. : 74469-00-4
Shelf Life : 3 years

Technical Specification

Appearance : White to off-white powder
Solubility : 33.3 mg soluble in 1 mL of water
pH : 8.00 - 10.00
Cultural response : Cultures conditions - Incubation period (5wks), Relative humidity (60±2%), Temperature (25±2°C), Photoperiod Day: Night in hours (16:8)
Shoot culture : No structural deformity observed.
Callus culture : No necrotic tissues, actively growing callus, no toxicity to callus.

Risk and Safety Information

Storage Temperature(°C) : Store at 2 - 8° C

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods