



Product Information

PR15

8-Strip PCR Tube Flat caps - Clear

Product Number

Packing

PR15 : 1 x 300 No.

Product Information

Product Code : PR15
Product Name : 8-Strip PCR Tube Flat caps - Clear

Technical Specification

Description : Autoclavable PCR tube flat caps connected together with in 8-strip format. It is made up of polypropylene material.
Clarity : Clear without any particles and fibers.
Height of Cap : 3.0 - 4.0 mm
Outer mouth Diameter of cap : 5.0 - 6.0 mm
Temperature tolerance test : Withstand steam sterilization at 121°C for 30 minutes and not to be used in dry heat sterilization.
Product alignment : 8- strip PCR tube flat caps connected to each other in straight line.

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein