

## PR19 96 Well Non-skirt, DNase RNase Free, Printed Alphanumeric PCR Blocks

Product Number	Packing
PR19	: 1 x 25 No.

### Product Information

---

Product Code	: PR19
Product Name	: 96 Well Non-skirt, DNase RNase Free, Printed Alphanumeric PCR Blocks

### Technical Specification

---

Description	: Neutral colour sterile non skirt PCR plate having alphanumerical printed marking.
Clarity	: Clear without any scratches and spots.
Well Format	: 96 well
Capacity	: 0.2 mL
DNase / RNase	: None Detected
Product alignment	: Proper alphanumeric marking and skirting.

### Transport Information

---

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods