

Disposable Serological Pipette (Narrow tip)

Product Code: PW1197

Product description:

PW1197 is a narrow tip, disposable serological pipette ideal for accurate liquid transfer or mixing. It has 50ml volume capacity.

Properties:

- Material: **Polystyrene or polyethylene**
- Volume capacity: **50ml**
- Graduations: **5/10**
- Colour code: **Purple**
- Sterile : **Yes**

Features:

- Made of high grade polystyrene or polyethylene for clear observation and minimal liquid attachment to the pipette surface
- Colour coded ring for easy identification
- Calibrated graduations for accurate dispensing within $\pm 2\%$
- Bidirectional graduations on the pipette provide added applicability
- Negative graduations allows additional working volume
- Supplied with a filter plug
- Leakage tested
- Universal design which makes the pipette suitable for most brands of pipette aid
- Sterilized by gamma irradiation
- Non-pyrogenic

Packaging:

- Individually packed in easy-to-open paper/plastic bag to assure sterility
- Every pack is properly labeled with lot number for quality traceability

Ordering information:

Code	Quantity/Pack	Total no. of Pack
PW1197-100X1no	100no	1