

## Disposable Petri Plates

**PW001, PW002, PW006, PW007, PW049, PW050, PW054, PW1128, PW1132, PW1132G, PW1134, PW1146, PW1151, PW1155, PW1173, PW1187, PW1190, PW1199, PW1199G, PW1223, PW1226, PW1232 & PW1244**

Transparent, disposable round plates with excellent optical clarity and strength. Petri plate recommended for laboratory use. For invitro use only.

**Application :** Bioassays ,Antimicrobial susceptibility tests , General lab purposes and Observation of plant germination.

| Product Name                                      | Product Code | Description  | Size (mm)      | Packing (Nos)  |
|---|--------------|--|----------------|--|
| <b>Sterile Disposable Petri Plate.</b>            | <b>PW001</b> | Optically clear, transparent disposable Petri plate  | 90X15          | 1X100 ,<br>1X540 &<br>20X540<br>Individually packed. |
| <b>Sterile Disposable Petri Plate.</b>            | <b>PW002</b> | Optically clear, transparent disposable Petri plate.   | 100X15         | 1X100 &<br>1X500<br>Individually packed.             |
| <b>Sterile Disposable Scored Top Petri Plate.</b> | <b>PW006</b> | Optically clear, transparent disposable and useful for culturing work.   | 90X15          | 1X100 &<br>1X500                                     |
| <b>Sterile Disposable Scored Top Petri Plate.</b> | <b>PW007</b> | Optically clear, transparent disposable and useful for culturing work.   | 65X12          | 1X100 &<br>1X500<br>Individually packed.             |
| <b>Sterile Disposable Petri Plate.</b>            | <b>PW049</b> | Optically clear, transparent disposable and useful for carrying out bioassay.                                    | 110X25         | 1X100 &<br>1X250                                     |
| <b>Sterile Disposable Square Petri Plate.</b>     | <b>PW050</b> | Optically clear, transparent disposable and useful for carrying out bioassay and antibacterial sensitivity test. | 120X120<br>X15 | 1X100  |

|  |                |  |         |  |
|--|----------------|--|---------|--|
| <b>Sterile Disposable Petri Plate. (Bi-Petri Plates)</b> | <b>PW054</b>   | Divided petri plates. Optically clear, transparent disposable having partition and useful for culturing work.    | 90X15   | 1X100 & 1X500                                    |
| <b>Sterile Disposable Scored Top Petri Plate.</b>        | <b>PW1128</b>  | Optically clear, transparent disposable petri plates made of polystyrene and useful for culturing work.          | 75X12   | 1X100  |
| <b>Sterile Disposable Petri Plate.</b>                   | <b>PW1132</b>  | Optically clear, transparent disposable.   | 90X15   | 10 X10, 54 X10 & 1080 X 10<br>Pack of 10 Plates. |
| <b>Sterile Disposable Petri Plates.</b>                  | <b>PW1132G</b> | Optically clear, transparent disposable and gamma irradiated.  | 90X15   | 10 X10, 54X10 & 200X10<br>Pack of 10 plates.     |
| <b>Sterile Disposable Petri Plate.</b>                   | <b>PW1134</b>  | Optically clear, transparent disposable petri-plates.  | 100X15  | 10X10 & 50X10<br>Pack of 10 plates.              |
| <b>Sterile Disposable Petri Plate.</b>                   | <b>PW1146</b>  | Optically clear, transparent disposable and useful for carrying out bioassay and antibacterial sensitivity test. | 120X15  | 1X100 & 1X250                                    |
| <b>Sterile Disposable Petri Plate.</b>                   | <b>PW1151</b>  | Optically clear, transparent disposable and useful for carrying out bioassay and antibacterial sensitivity test. | 110X15  | 1X100 & 1X250                                    |
| <b>Sterile Disposable Petri Plate.</b>                   | <b>PW1155</b>  | Optically clear, transparent disposable petri plate.   | 55X12   | 5X25 & 10X25.<br>Pack of 25 plates.              |
| <b>Sterile Disposable Petri Plate, Triple Pack.</b>      | <b>PW1173</b>  | Optically clear, transparent disposable petri plate.   | 90X15   | 10X10 & 50X10                                    |
| <b>Sterile Disposable Scored Top Petri Plate.</b>        | <b>PW1187</b>  | Optically clear, transparent disposable and useful for culturing work.   | 55 X 12 | 1X100 & 1X500                                    |
| <b>Sterile Disposable Petri Plate.</b>                   | <b>PW1190</b>  | Optically clear, transparent disposable petri plate.   | 90X25   | 1X100 & 1X250<br>Individually                    |

|   |                |  |        |  |
|---|----------------|--|--------|--|
|   |                |  |        | packed.                                      |
| <b>Sterile Disposable Petri Plate.</b>              | <b>PW1199</b>  | Optically clear, transparent disposable petri plate.   | 90X15  | 5 X20, 25X20 & 500X20.<br>Pack of 20 plates. |
| <b>Sterile Disposable Petri Plate.</b>              | <b>PW1199G</b> | Optically clear, transparent disposable and gamma irradiated petri plates.   | 90X15  | 5 X20 , 25X20 & 500X20<br>Pack of 20 plates. |
| <b>Sterile Disposable Petri Plate.</b>              | <b>PW1223</b>  | Optically clear having 6 vent.   | 55X14  | 5X15 & 10X15<br>Pack of 15 plates.           |
| <b>Sterile Disposable Petri Plate.</b>              | <b>PW1226</b>  | Optically clear, transparent disposable petri plates. Useful for carrying out bioassay and antibacterial sensitivity test. | 150X20 | 1X100  |
| <b>Sterile Disposable Scored Top Petri Plate.</b>   | <b>PW1232</b>  | Optically clear, transparent disposable and useful for culturing work.   | 65X12  | 10X10 & 50X10<br>Pack of 10 plates.          |
| <b>Sterile Disposable Petri Plate.(Polystyrene)</b> | <b>PW1244</b>  | Optically clear, transparent disposable petri plates made of polystyrene.  | 110X40 | 1X100 & 1X250<br>Individually packed.        |

**Warning :** Avoid dry heat sterilization.

**Disposal :** User must ensure safe disposal by autoclaving at a temperature of 121°C (15 lbs) for 15 minutes of used medium (bacteriological materials e.g. culture media with pathogens or spores) in petri plates / unused preparations of this product on completion of work to avoid any contagion.

**Product Features :**

- Optically Clear for Easy Reading.
- Packaged sterile.
- These Disposable Polystyrene Petri Dishes are supplied in a variety of sizes.

**Disclaimer :**



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.