

**Thumb Press Dispensing Dropper / pasture Pipette Plastic**      **PW136, PW1144, PW1145, PW1201, PW1206, PW1212**

**Thumb Press Dispensing Dropper / pasture Pipette** : Plastic Droppers /pipettes made from non-toxic polyethylene plastic. They are semi transparent and have gradations that continue down into the tips. Graduation ranges from 1.00 ml and can easily dispense 1 drop at a time.

These dispensing droppers /pipettes have a top where the user holds his thumb over the top to create the vacuum and seal the liquid in the Dropper / pipette. Blow-out dispensing droppers / pipettes, allow you to blow into this open end to get the last bits of liquid in the droppers / pipette into your receiving container, for greater accuracy.

Blow-out Dispensing droppers / pipettes are clearly marked with a frosted thin rings around the neck.

**How It works :**

- Hold the dispensing dropper / pipette with Pipette Tips vertically at eye level after you have drawn up your liquid.
- Find the level of the liquid inside the dispensing dropper / pipette. The surface of the liquid will look as though it curves. This is called a “meniscus curve”.
- Find the measurement line that meets the bottom of the meniscus curve.
- Count up the number of lines from the nearest whole number below the bottom of your meniscus curve. Record your reading with the appropriate units.

**Application:** Use in plant tissue culture work, Molecular Biology work, Laboratories, Universities and other Industries .

Product Name	Product Code	Description	Capacity (ml)	Material
<b>Thumb Press Dispensing Dropper / pasture Pipette</b>	<b>PW136</b>	Thumb Press Dispensing Dropper / pasture Pipette	3.00	Polyethylene
	<b>PW1144</b>	Thumb Press Dispensing Dropper / pasture Pipette	1.00	Polyethylene
	<b>PW1145</b>	Thumb Press Dispensing Dropper / pasture Pipette	2.00	Polyethylene
	<b>PW1201</b>	Thumb Press Dispensing Dropper / pasture Pipette having drop volume 25.00 µl	1.00	Polyethylene
	<b>PW1206</b>	Thumb Press Dispensing Dropper / pasture Pipette having drop volume 20.0 µl	2.00	Polyethylene
	<b>PW1212</b>	Thumb Press Dispensing Dropper / pasture Pipette And they are available in individual sterile pack having drop volume 3.00 µl.	3.00	Polyethylene

**Product Features :**

- They are graduated pipettes.
- Non-Toxic.
- Available in natural Colour.
- Pack size : PW136, PW1144 & PW1145 -1 X50 Nos and  
PW1201 , PW1206 & PW1212 - 1 X 50 Nos. & 1 X 500 Nos.

**Disclaimer :**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.