



# Product Information

**PW1215**

**100 wells polypropylene cryogenic storage box**

---

**Product Number****Packing**

PW1215 : 1 X 5 No.

---

**Product Information**

Product Code : PW1215  
Product Name : 100 wells polypropylene cryogenic storage box

---

**Technical Specification**

Description : Colourless plastic, autoclavable cryogenic storage box with alphanumerically labeled places for eppendorf tubes or cryo tubes, made up of polypropylene material, having translucent lid.

Dimension : 55 x 141 x 141 mm

Clarity : Clear translucent without any spots and breakage.

Capacity : 100 wells

Temperature tolerance test : Withstands steam sterilization at 121°C for 30 minutes and not to be used in dry heat sterilization.

---

**Transport Information**

ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein