



Product Information

PW1255

0.2 mL PCR Tubes, Thin walled

Product Number**Packing**

PW1255 : 1 X 1000 No

Product Information

Product Code : PW1255
Product Name : 0.2 mL PCR Tubes, Thin walled

Technical Specification

Description : Autoclavable PCR tubes having conical base with flat lid. It is DNA/RNase free tubes and made up of polypropylene material.
Clarity : Clear without any particles and fibers.
Height of Tube : 19.5 - 20.5 mm
Height of Cap : 3.5 - 4.5 mm
Inner mouth Diameter of tube : 5.0 - 6.0 mm
Outer mouth Diameter of tube : 6.5 - 7.5 mm
Outer mouth Diameter of cap : 5.0 - 6.0 mm
Capacity : 0.2 mL
Temperature tolerance test : Withstand steam sterilization at 121°C for 30 minutes and not to be used in dry heat sterilization.
Product alignment : Lid closer of the tube is tight.

Transport Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein