

## Serum Vial Sample Storage Box and Stand

PW145

Storage box consists of a base and insert. The insert is removable for storage versatility. The storage box stores 96 serum vials. The box and lid are one way on/off fit preventing orientation errors. It stores Ninety six - 3.0 ml vials capacity in 12 X 8 format.

The rack enables single-handed vial manipulation. The top of each opening interlocks with the neck of each vial permitting closures to be removed or tightened with one hand. Racks are stackable when empty or full.

**Application:** Microbiology laboratory, Molecular Biological laboratory, Hospitals, Cell Culture and Various other laboratories.

| Product Name                             | Product Code | Description   |
|--|--------------|---|
| Serum Vial Sample Storage Box and Stand. | PW145        | 96 vials fitted in the box. Disposable transparent screw capped vials with 3.0 ml capacity. |

### Product Features :

- For use of storage vials.
- Stores 96 vials in 12 X 8 Format.
- Stackable.
- Pack size : 1 X 96/ pack.

### Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.