

Coplin Jar

PW354

They are Autoclavable, Air tight, water tight and Dome style screw capped with self standing capability. Three sides grooved bottom and two sides to protect glass slides, wider mouth provided to put the slides properly and draw the slides without using any forceps, unbreakable, non reactive with any solvent , available in white, grey, yellow and red colours.

Application : Microbiology laboratory, Plant tissue culture work, Hospitals, Pathology and various other laboratories.

Product Name	Product Code	Description	Materials
Autoclavable Coplin Jar	PW354	Air Tight, No moisture absorption, Unbreakable, Self standing jar to protect glass slides.	Polypropylene Co-Polymer (PPCP)

Product Features :

- Autoclavable.
- Unbreakable
- Non reactive with solvents.
- Pack Size : 100 Nos/Pack.

Disclaimer :



User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.