

RM10793 Poly(ethylene glycol) MW 10000

Product Number Packing

RM10793 : 500G

Product Information

Product Code : RM10793
Product Name : Poly(ethylene glycol) MW 10000
Synonym : PEG
Molecular Formula : $H(OCH_2CH_2)_nOH$
Average Molecular Weight : ~10,000
CAS No. : 25322-68-3
HS Code : 3907 20 11
Shelf Life : 4 years

Technical Specification

Appearance : White crystals or powder or flakes
Solubility : 33.3 mg soluble in 1 mL of water
pH (5% in water at 25°C) : 5.00 - 7.00
Melting range : 58 - 64°C

Risk and Safety Information

WGK : 1
Storage Temperature(°C) : Store below 30°C

Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods