

RM10799 Poly(3-hexylthiophene-2,5-diyl)**Product Number Packing**

RM10799 : 250MG

Product Information

Product Code : RM10799
Product Name : Poly(3-hexylthiophene-2,5-diyl)
Molecular Formula : $(C_{10}H_{14}S)_n$
CAS No. : 156074-98-5
HS Code : 3212 90 00
Shelf Life : 4 years

Technical Specification


Appearance : Brown to black powder or solid or fibre
Solubility : 33.3 mg soluble in 1 mL of chloroform
Melting range : 233 - 243°C
Assay (Based on Trace Metal analysis): min. 99.99%

GHS Safety Information

Hazard Statement(s) : H315-H319-H335
Precautionary Statement(s) : P261-P305+P351+P338
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 36/37/38
S-Phrase(s) : 26
WGK : 3
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods