

RM10844 Carboxymethylcellulose sodium salt, Medium viscosity (CMC)

Product Number Packing

RM10844 : 500G

Product Information

Product Code : RM10844
Product Name : Carboxymethylcellulose sodium salt, Medium viscosity (CMC)
Synonym : Sodium carboxymethylcellulose
CAS No. : 9004-32-4
HS Code : 3912 31 00
Shelf Life : 4 years

Technical Specification

Appearance : White to yellow to light tan powder
Solubility : 10 mg soluble in 1 mL of water
pH (1% in water at 25°C) : 6.50 - 8.00
FTIR : Matches with the standard pattern
Loss on drying (at 110°C, 2 hrs.) : <= 15.00%
Viscosity (2% in water at 25°C) : 400 - 800 cps

Risk and Safety Information

WGK : 1
Storage Temperature(°C) : Store below 30°C

Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods