

RM10905 Lipase, From Aspergillus niger

Product Number Packing

RM10905 : 100G
RM10905 : 25G

Product Information

Product Code : RM10905
Product Name : Lipase, From Aspergillus niger
CAS No. : 9001-62-1
EC No. : 232-619-9
HS Code : 3507 90 90
Enzyme Classification No. : 3.1.1.3
Shelf Life : 3 years

Technical Specification

Appearance : Off-white lyophilized powder
Solubility : 2 mg soluble in 1 mL of water forming hazy faintly yellow coloured solution with insoluble particles
Activity : ~ 50 U/mg
Unit definition : 1U corresponds to the amount of enzyme which hydrolyzes 1µmol acetic acid per minute at pH 7.4 and 40°C

Risk and Safety Information

WGK : 1
RTECS : TO9776500
Storage Temperature(°C) : Store at 2 - 8°C

Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods