

RM1734 Octyl sulphate sodium salt, For HPLC**Product Number Packing**

RM1734	: 1G
RM1734	: 5G


Product Information

Product Code	: RM1734
Product Name	: Octyl sulphate sodium salt, For HPLC
Synonym	: Sodium octyl sulphate, For HPLC
Molecular Formula	: C ₈ H ₁₇ NaO ₄ S
Molecular Weight	: 232.27
CAS No.	: 142-31-4
EC No.	: 205-535-5
HS Code	: 2920 90 10
Shelf Life	: 4 years


Technical Specification

Appearance	: White to off-white hygroscopic powder or solid
Solubility	: 33.3 mg soluble in 1 mL of water
Melting range	: 195 - 210°C (dec.)
UV Absorption 1	: ≤ 0.15 (at 200 nm, 0.005M/L in water)
UV Absorption 2	: ≤ 0.1 (at 210 nm, 0.005M/L in water)
UV Absorption 3	: ≤ 0.06 (at 220 nm, 0.005M/L in water)
UV Absorption 4	: ≤ 0.04 (at 230 nm, 0.005M/L in water)
UV Absorption 5	: ≤ 0.02 ((at 260 nm, 0.005M/L in water)
Assay (Hyamine Titration)	: min. 99.00 %

GHS Safety Information

Hazard Statement(s)	: H315-H319-H335
Precautionary Statement(s)	: P261-P305+P351+P338
Signal Word	: Warning
Hazard Pictogram(s)	

Risk and Safety Information

R-Phrase(s)	: 36/37/38
S-Phrase(s)	: 26-36
WGK	: 3
RTECS	: WT 1175000
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	 Irritant



Product Information

Transport Information

Marine Pollutant	:	No
ADR/RID	:	Not Dangerous Goods
IMDG	:	Not Dangerous Goods
IATA	:	Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein