



# Product Information

**RM2662 Biogel P-10, 100-200 mesh**

---

<b>Product Number</b>	<b>Packing</b>
-----------------------	----------------

RM2662	: 100G
RM2662	: 25G

---

**Product Information**

Product Code	: RM2662
Product Name	: Biogel P-10, 100-200 mesh
Molecular Weight	: 1500 - 20000
Other Information	: Polyacrylamide beads
Shelf Life	: 4 years

---

**Technical Specification**

Appearance	: White powder.
Solubility	: Insoluble in water.
Particle size	: 100 - 200 mesh

---

**Risk and Safety Information**

Storage Temperature(°C)	: Store below 30°C
-------------------------	--------------------

---

**Transport Information**

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein