



# Product Information

**RM3357 DEAE Dextran**

---

**Product Number Packing**

RM3357 : 100G  
RM3357 : 25G

---

**Product Information**

Product Code : RM3357  
Product Name : DEAE Dextran  
Average Molecular Weight : 5,00,000  
CAS No. : 9015-73-0  
Shelf Life : 4 years

---

**Technical Specification**

Appearance : White to cream hygroscopic powder  
Solubility : 100 mg soluble in 1 mL of water  
pH (5% in water at 25°C) : 4.00 - 7.00  
Loss on drying (at 105°C, 2hr) : <= 7.00%  
Nitrogen content : 3.00 - 3.50%

---

**Risk and Safety Information**

WGK : 3  
Storage Temperature(°C) : Store below 30°C

---

**Transport Information**

Marine Pollutant : -  
ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

-----  
The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein