

**RM3882 Ammonium acetate, ACS**

**Product Number Packing**

---

RM3882 : 500G  
RM3882 : 5KG

**Product Information**

---

Product Code : RM3882  
Product Name : Ammonium acetate, ACS  
Molecular Formula :  $C_2 H_7 N O_2$   
Molecular Weight : 77.08  
CAS No. : 631-61-8  
EC No. : 211-162-9  
HS Code : 2915 29 00  
Shelf Life : 4 years

**Technical Specification**

---

Appearance : Colorless or white deliquescent crystals.  
Solubility : 33.3 mg soluble in 1 mL of water  
pH (5% in water at 25°C ) : 6.70 - 7.30  
Chloride (Cl) :  $\leq 0.0005\%$   
Heavy metals (as Pb) :  $\leq 0.0005\%$   
Iron (Fe) :  $\leq 0.0005\%$   
Sulphate ( $SO_4$ ) :  $\leq 0.001\%$   
Nitrate ( $NO_3$ ) :  $\leq 0.001\%$   
Insoluble matter :  $\leq 0.005\%$   
Residue after ignition :  $\leq 0.01\%$   
Assay (NaOH back Titration) : min. 97.00 %

**GHS Safety Information**

---

Hazard Statement(s) : H315-H319-H335  
Precautionary Statement(s) : P261-P305+P351+P338  
Signal Word : Warning  
Hazard Pictogram(s)



**Risk and Safety Information**

---

RTECS : AF3675000  
Storage Temperature(°C) : Store below 30°C

**Transport Information**

---



# Product Information

## Transport Information

---

Marine Pollutant	:	No
ADR/RID	:	Not Dangerous Goods
IMDG	:	Not Dangerous Goods
IATA	:	Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein