



Product Information

RM8036 Carboxy vinyl polymer ETD-2001

Product Number	Packing
----------------	---------

RM8036	: 100G
--------	--------

Product Information

Product Code	: RM8036
Product Name	: Carboxy vinyl polymer ETD-2001
Synonym	: Carbopol [®] ETD 2001
Other Information	: 1. Crosslinked polyacrylate polymer. 2. Used as thickening/gelling agent in personal care application. 3. Key : ® : Registered trademark of The Lubrizol Corporation
Shelf Life	: 4 years

Technical Specification

Appearance	: White powder.
Solubility	: Reacts with fatty amines to form thick and stable emulsions of oils in water

Risk and Safety Information

Storage Temperature(°C)	: Store below 30°C
-------------------------	--------------------

Transport Information

Marine Pollutant	: -
ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein