



Product Information

RM8164 **Dextran, MW 50,00,000 - 400,00,000**

Product Number	Packing
----------------	---------

RM8164	: 25G
--------	-------

Product Information

Product Code	: RM8164
Product Name	: Dextran, MW 50,00,000 - 400,00,000
Molecular Formula	: $(C_6H_{10}O_5)_n$
Molecular Weight	: 50,00,000 - 400,00,000
CAS No.	: 9004-54-0
EC No.	: 232-677-5
HS Code	: 3913 90 00
Shelf Life	: 4 years

Technical Specification

Appearance	: White to off-white powder or flakes
Solubility	: 10 mg soluble in 1 mL of water giving hazy solution
FTIR (KBr disc)	: Matches with the standard pattern

Risk and Safety Information

WGK	: 2
Storage Temperature(°C)	: Store below 30°C

Transport Information

Marine Pollutant	: No
ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein