

**RM857                      Dextran sulfate sodium salt, From Leuconostoc spp.****Product Number****Packing**

RM857	: 25G
RM857	: 5G

**Product Information**

Product Code	: RM857
Product Name	: Dextran sulfate sodium salt, From Leuconostoc spp.
Molecular Weight	: 5,00,000
CAS No.	: 9011-18-1
HS Code	: 3913 90 00
Shelf Life	: 4 years

**Technical Specification**

Appearance	: White to light yellow powder or solid
Solubility	: 100 mg soluble in 1 mL of water
pH (2% in water at 25°C)	: 6.0 - 8.0
Loss on drying (at 105°C, 2 hr)	: ≤ 10.00%
Specific rotation	: +75.00° to +105.00° (c = 4% in water at 20°C)

**Risk and Safety Information**

WGK	: 2
RTECS	: HH9290000
Storage Temperature(°C)	: Store at 2 - 8°C

**Transport Information**

Marine Pollutant	: No
ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods