



## Sterile Discs\*

SD067

Sterile discs can be used for testing antibiotic susceptibility testing by loading antibiotic of choice.

### Appearance:

Sterile filter paper discs of 6mm diameter

### Cultural response:

Sterility of the discs checked for 14 days in Alternative Thioglycollate Medium (M010) at 30-35°C and Soyabean Casein Digest Medium (M011) at 20-25°C.

Sample	Results in Alternative Thioglycollate Medium (M010)	Results in Soyabean Casein Digest Medium (M011)
Test	Clear, No growth observed	Clear, No growth observed
Negative Control	Clear, No growth observed	Clear, No growth observed
Positive control	Luxuriant growth observed	Luxuriant growth observed

### Storage and Shelf-life:

Store at 2°C to 8°C. Use before expiry date on the label.

\* Not for Medicinal Use