

---

## Pyridoxal hydrochloride

**TC164**

### Cell Culture Tested

#### Product Identifier

---

CAS No.	:	65-22-5
EC No.	:	200-602-5
Molecular Formula	:	C <sub>8</sub> H <sub>9</sub> NO <sub>3</sub> .HCl
Molecular Weight	:	203.62
HS Code	:	2936 25 00
Storage	:	2 - 8°C
Shelf life	:	3 years

#### Technical Specification

---

Appearance	:	White to pale brown crystals or powder or solid
Solubility	:	Clear solution at 100 mg soluble in 1 mL of water
pH (1% in water at 25°C)	:	2.00 - 4.00
FTIR	:	Matches with the standard pattern
Assay (NT/AT)	:	min. 99.00%

#### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	UV1225000
WGK	:	2