

TC303 Dispase®, Cell Culture Tested**Product Number Packing**

TC303 : 100MG

Product Information

Product Code : TC303
Product Name : Dispase®, Cell Culture Tested
CAS No. : 42613-33-2
EC No. : 255-914-4
HS Code : 3507 90 90
EC Index No. : 647-014-00-9
Other Information : Key : ® : Dispase is a registered trademark of Godi Shusei Co. Ltd., Tokyo, Japan
Shelf Life : 3 years

Technical Specification

Appearance : White to brown powder.
Cell Culture Test : Passes
Enzymatic Activity : min. 0.50 units/mg

GHS Safety Information

Hazard Statement(s) : H315-H319-H334-H335
Precautionary Statement(s) : P261-P305+P351+P338-P342+P311
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 36/37/38-42
S-Phrase(s) : 22-24-26-36/37
WGK : 3
Storage Temperature(°C) : Store at 2 - 8°C
Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods