
17-AAG, Cell Culture Tested

TC409

Product Identifier

CAS No.	:	75747-14-7
Molecular Formula	:	$C_{31}H_{43}N_3O_8$
Molecular Weight	:	585.69
Synonym	:	17-(Allylamino)-17-demethoxygeldanamycin
HS Code	:	2941 90 90
Storage	:	2 - 8°C
Shelf life	:	TC409

Technical Specification

Appearance	:	Purplish to purplish red solid
Solubility	:	33.3 mg soluble in 1 mL of methanol
Assay (HPLC)	:	98.00 - 102.00%
Cell Culture Test	:	Passes

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3