

## Dulbecco's Phosphate Buffered Saline 1X

With 1000mg/L Glucose, 36mg Sodium pyruvate per litre, Calcium chloride, Magnesium sulphate, 50mg Streptomycin sulphate per litre, 100mg Kanamycin sulphate per litre Without Phenol red

**Product Code: TL1123**

### Product Description :

All media used in tissue culture have a basis of a synthetic mixture of inorganic salts known as a physiological or balanced salt solution (BSS). All the physiological salt solutions have been derived from the salt solution originally described by Sydney Ringer (1885). The first balanced salt solution to be developed specifically for supporting the metabolism of mammalian cells was Tyrode's solution. Since then many modifications have been done to obtain better buffering salt solutions and to prevent calcium precipitation.

The function of a salt solution is:

- To maintain the medium within physiological pH range.
- To maintain intracellular and extra cellular osmotic balance.
- Modified with a carbohydrate, such as glucose serves as an energy source for cell metabolism.

TL1123 is 1X Dulbecco's Phosphate Buffered Saline with 1000mg per litre D-glucose, 36mg per litre sodium pyruvate, calcium chloride, magnesium sulphate, 50mg per litre streptomycin sulphate, 100mg per litre kanamycin sulphate. It does not contain phenol red. It is frequently used with primary cells, where contamination risk is high.

### Composition :

Ingredients	mg/L
<b>INORGANIC SALTS</b>	
Calcium chloride dihydrate	133.000
Disodium hydrogen phosphate	1150.000
Magnesium chloride hexahydrate	100.000
Potassium chloride	200.000
Potassium phosphate, monobasic	200.000

Sodium chloride	8000.000
<b>OTHERS</b>	
D - Glucose	1000.000
Kanamycin sulphate	100.000
Sodium pyruvate	36.000
Streptomycin sulphate	50.000

### Quality Control:

#### Appearance

Colorless, clear solution

#### pH

6.90 -7.50

#### Osmolality in mOsm/Kg H<sub>2</sub>O

300.00 -340.00

#### Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

#### Toxicity test

Passes

#### Endotoxin content

NMT 1EU/ml

### Storage and Shelf Life:

Store at 2-8°C away from bright light.

Shelf life is 24 months.

Use before expiry date given on the product label.

Revision: 02 / 2019

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.